

## DISCOVERY PICO LASER TATOO REMOVAL PRE/POST TREATMENT

1 of 2

## BEFORE DISCOVERY PICO LASER TATTOO REMOVAL:

## TATOO REMOVAL:

- TREATMENTS: The number of treatments depends on each individual's unique circumstances, but the Discovery Pico achieves optimum outcome in fewer sessions than almost any other tattoo removal machine.
- REDUCED PAIN: The Discovery Pico is designed to reduce the discomfort and recovery time often associated with laser tattoo removal. Skin heals quickly and evenly, and there is little risk of permanent scarring or discoloration of the skin.
- HIGHLY EFFECTIVE: The Discovery Pico laser outclasses conventional picosecond lasers by using the highest power possible, 3 wavelengths, and 4 emission modes to fragment even the deepest pigments precisely and evenly. It even works on hard-to-treat tattoo ink colors and pigments like light blue and green.

## FREQUENTLY ASKED QUESTIONS ABOUT LASER TREATMENTS:

- Is a Quanta Discovery PICO treatment right for me? Quanta Discovery PICO lasers are designed to deliver a wide range of safe and effective treatments for all skin types. These innovative lasers can minimize benign pigmented lesions, remove acne scarring and unwanted tattoos, and smooth fine lines and wrinkles by replenishing collagen production
- Will my Discovery PICO treatment hurt? The Discovery PICO lasers are designed to deliver comfortable treatments, while reducing recovery time. The lasers treat pigmentation as evenly and precisely as possible to eliminate hot spots and promote natural healing.
- How many treatment sessions will I need? The number of treatments needed varies depending on each patient's desired outcome. During your consultation, we will be able to determine your treatment plan.
- What Sets the Discovery PICO Series Apart: Minimal discomfort, Speedy recovery time, Unparalleled precision & Renewal of confidence.

Discovery PICO represents the second generation of picosecond lasers that due to its proprietary technology, is the most innovative medical laser available today in the market. It has an Nd:YAG double-wavelength laser source emitting at both 1064nm and 532 nm.

Discovery Pico Plus is the first picosecond yag laser with the same specifications of Discovery Pico combined with a ruby laser able to delivers up to 1.2J of energy and not just few hundred mJ such as other "red" lasers. The 694 nm of Discovery Pico Plus are generated by the most advanced and powerful ruby laser available in the market: very high energy per pulse in Q-Switched or Photo Thermal pulses, big spot sizes and fast repetition rates make this laser, together with the 532 and 1064 nm which work also in picosecond regime.

Discovery PICO series can treat different tattoo colors, superficial and deep pigmented lesions by using Q-Switched and Photo-Thermal pulse emission modes, as an exclusive feature of this laser, with possibilities to treat aged skin effectively.

The innovation introduced by Discovery PICO series is the level of peak power through ultrashort pulses that is twice more powerful than any other Picosecond laser for medical applications. This has enhanced the treatment results increasing the efficacy on hard-to-treat tattoo inks with a general reduction of the delivered average energy dose compared to traditional Q-Switched lasers.

Discovery PICO series is indicated for the use of:

- Tattoo Removal
- Benign Pigmented lesions Removal
- Hair Lightening
- Skin Rejuvenation
- Foot Skin and Nail condition